

Personal inductive neck loop

Inductive plastic neck loop with 3.5mm jack plug


The inductive neck loop can be used with any equipment with a suitable headphone socket eg. personal radios, minidisks, CD players, cassettes, talking book machines etc.

Simply plug directly into the headphone socket and switch the hearing aid to the 'T' setting.

The loop works equally well when used with mono or stereo equipment and can be worn under clothes without loss of performance.

With 3.5mm stereo plug & 6.3mm stereo adaptor.

*Can also be used with 3 metre stereo extension lead
(Pt No. ★A121A★)*



Plastic neck loop with 3.5mm plug
(Pt No. 95THINSB)



CONNEVANS LIMITED

54 Albert Road North, Reigate, Surrey RH2 9YR

Helpline 01737 247571

help@connevens.com