

Australia

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,400	15	H02
173,640	26	H06
173,990	69	H33
174,120	41	H16
174,415	51	H20
174,600	55	H06
174,770	61	H08
175,700	2	H27
176,000	11	H30

Austria

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
174,030	73	
174,070	39	H15
174,120	41	H16
174,185	43	H17
174,270	46	H18
174,360	49	H19
174,415	51	H20
174,600	55	H06
174,640	56	
174,680	74	
174,720	75	
174,800	62	
174,880	70	
174,920	A0	
175,000	67	

Belgium

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
181,900	67	H24
183,825	28	
183,875	29	
184,180	68	

Canada

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
216,01250	1	N01
216,03750	2	N02
216,08750	4	N04
216,11250	5	N05
216,18750	8	N08
216,21250	9	N09
216,28750	12	N12
216,31250	13	N13
216,38750	16	N16
216,41250	17	N17
216,43750	18	N18
216,51250	21	N61
216,53750	22	N62
216,58750	24	N64
216,61250	25	N65
216,68750	28	N68
216,71250	29	N69
216,78750	32	N72
216,81250	33	N73
216,88750	36	N76
216,91250	37	N77
216,96250	39	N79
216,98750	40	N80
216,07500	42	
216,22500	45	
216,37500	48	
216,52500	51	
216,57500	52	N52
216,67500	54	
216,72500	55	N55
216,82500	57	N57
216,97500	60	N60
169,4450	B1	
169,5050	B2	
170,2450	B3	
170,3050	B4	
171,0450	B5	
171,1050	B6	
171,8450	B7	
171,9050	B8	

Croatia

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,4	15	H02
173,465	19	H03
173,545	23	H04
173,640	26	H06
173,695	28	H11
173,755	31	
173,950	36	H14
174,120	41	H16
174,415	51	H20
174,600	55	H06
174,770	61	H08

Czech Republic

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
174,120	41	H16
174,175	77	H77
174,275	87	H87
174,325	78	H78
174,425	89	H89
174,600	55	H06
175,000	67	H67

Denmark

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
169,675	P6	H45
169,725	P7	H99
169,775	P8	H46
169,825	P9	
173,400	15	H02
173,465	19	H03
173,640	26	H05
173,695	28	H11
173,813	91	H35
173,825	33	H13
173,963	93	H34
180,500	70	
187,500	71	J20
208,500	61	
215,500	62	

Finland

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,990	69	H33
174,175	77	H77
174,275	87	H87
174,325	78	H78
174,425	89	H89
174,475	79	H79
174,600	55	H06
175,000	67	H67
175,700	2	H27
175,850	97	
175,950	10	
176,000	11	H30

France

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
175,500	1	H25
175,600	95	H26
175,700	2	H27
175,800	9	H28
175,900		H29
176,000	11	H30
183,925	30	
183,975	31	
184,125	32	
184,175	33	
184,510	69	

Germany

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,990	69	H33
174,075	84	H84
174,125	85	H85
174,175	77	H77
174,225	86	H86
174,275	87	H87
174,325	78	H78
174,375	88	H88
174,425	89	H89
174,475	79	H79
174,525	90	H90
174,600	55	H06
174,650	57	H57
174,700	59	H59
174,775	81	H81
175,000	67	H67

Italy

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,988	83	H96
173,990	69	H33
174,000	37	

Korea

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
217,425	69	
217,525	71	
217,675	74	
217,775	76	
217,925	79	

The Netherlands

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,4	15	H02
173,465	19	H03
173,545	23	H04
173,640	26	H06
173,695	28	H11
173,755	31	
173,950	36	H14
174,120	41	H16
174,415	51	H20
174,600	55	H06
174,770	61	H08

New Zealand

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,000	A1	H96
173,010	A8	H91
173,050	A2	
173,100	A3	H92
173,150	A4	
173,200	A5	H93
173,250	A6	
173,300	A7	H94
173,340	12	
173,350	13	H01
173,400	15	H02
173,430	16	
173,450	17	
173,460	18	
173,465	19	H03
173,490	20	
173,500	21	
173,520	22	
173,545	23	H04
173,550	24	
173,600	25	
173,640	26	H05
173,650	27	
173,695	28	H11
173,700	29	
173,750	30	
173,755	31	H95
173,775	72	H12
173,800	32	
173,813	91	H35
173,825	33	H13
173,850	34	
173,900	35	
173,950	36	H14
173,963	93	H34
173,990	69	H33
174,000	37	
174,050	38	
174,100	40	

Norway

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,8125	91	H35
173,8375	7	H91
173,8625	94	H92
173,8875	76	H93
173,9125	92	H94
173,9375	8	H95
173,9625	93	H34
173,9875	83	H96

Poland

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
174,120	41	H16
174,175	77	H77
174,275	87	H87
174,325	78	H78
174,425	89	H89
174,600	55	H06
174,850	63	
175,000	67	H67
175,700	2	H27
175,850	97	
175,950	10	
176,000	11	

South Africa

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,755	31	
173,950	36	H14
174,325	78	H78
174,425	89	H89
174,475	79	H79
174,770	61	H08

Spain

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
174,10000	40	H95
174,30000	47	
175,50000	1	H25
176,30000	4	
180,00000	13	

Sweden

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
169,6750	P6	H45
169,7250	P7	H99
169,7750	P8	H46
169,8250	P9	
173,4000	15	H02
173,8125	91	H35
173,9625	93	H34
173,9900	69	H33
174,1200	41	H16
174,1750	77	H77
174,2750	87	H87
174,3250	78	H78
174,4150	51	H20
174,6000	55	H06
174,7700	61	H08
174,8500	63	
175,0000	67	H67
175,7000	2	H27
175,8500	97	
175,9500	10	
176,0000	11	H30
180,4500	72	H31
181,2100	64	H21
187,4500	73	J19
208,0000	90	N23
216,0125	1	N01
216,2125	9	N09
216,4375	18	N18

Switzerland

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
174,120	41	H16
174,175	77	H77
174,275	87	H87
174,425	89	H89
174,600	55	H06
174,770	61	H08
175,000	67	H67
175,700	2	H27

Taiwan

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
216,0125	1	N01
216,2875	12	N12
216,3125	13	N13
216,4375	18	N18
216,5875	24	N64
216,7875	32	N72
216,9125	37	N77
216,9875	40	N80
216,075	42	
216,225	45	
216,375	48	
216,525	51	
216,675	54	
216,825	57	N57
216,975	60	N60

Turkey

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,990	69	H33
174,175	77	H77
174,275	87	H87
174,325	78	H78
174,425	89	H89
174,475	79	H79
174,600	55	H06
175,000	67	H67
175,700	2	H27
175,850	97	
175,950	10	
176,000	11	H30

United Kingdom

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
173,350	13	
173,400	15	H02
173,450	17	
173,500	21	
173,600	25	
173,650	27	
173,700	29	
173,750	30	
173,800	32	
173,850	34	
173,900	35	
173,950	36	
174,000	37	
174,050	38	
174,100	40	
174,150	42	
174,200	44	
174,250	45	
174,300	47	
174,350	48	
174,400	50	
174,450	52	
174,500	53	
174,550	54	
174,600	55	H06
174,650	57	
174,700	59	
174,750	60	
174,800	62	
174,850	63	
174,900	65	
174,950	66	
175,000	67	
175,050	71	

USA

List of total accessible channels		
Frequency	Oticon Channel	Phonak Channel
216,01250	1	N01
216,03750	2	N02
216,08750	4	N04
216,11250	5	N05
216,18750	8	N08
216,21250	9	N09
216,28750	12	N12
216,31250	13	N13
216,38750	16	N16
216,41250	17	N17
216,43750	18	N18
216,51250	21	N61
216,53750	22	N62
216,58750	24	N64
216,61250	25	N65
216,68750	28	N68
216,71250	29	N69
216,78750	32	N72
216,81250	33	N73
216,88750	36	N76
216,91250	37	N77
216,96250	39	N79
216,98750	40	N80
216,07500	42	
216,22500	45	
216,37500	48	
216,52500	51	
216,57500	52	N52
216,67500	54	
216,72500	55	N55
216,82500	57	N57
216,97500	60	N60
169,4450	B1	
169,5050	B2	
170,2450	B3	
170,3050	B4	
171,0450	B5	
171,1050	B6	
171,8450	B7	
171,9050	B8	