

Connevens Guide to Hearing Aid Direct Input Shoes – 1

Some hearing aids have a direct connection facility, usually referred to as direct input. This facility allows the direct connection of the hearing aid to other audio equipment via a connector usually referred to as a shoe; thus the name 'Hearing aid direct input shoe'. Connevens is able to provide, usually from stock, a wide range of hearing aid direct input shoes; many customers find the number of options overwhelming!

You can put your shoe against the page and identify most shoes without difficulty, making re-ordering easier for you and us!


If you require further help, please do phone our radio aid helpline for assistance when ordering. Please note that although we have usually shown two shoes in the picture to give different views, the shoes are all sold as singles.

VAT All direct input shoes: VAT relief is permitted subject to the receipt of a valid VAT declaration form.



Information on compatibility with ear level receivers

By each shoe we indicate whether it is compatible for use with 'ear level receivers' such as MicroMLxS, MLx, MLxS, Lexis, Amigo, Scola and fmGenie WDI units. For use with an MLxS, we also show whether the MLxS needs to be set to a crosswise or lengthwise orientation. For use with fmGenie WDI and MicroMLxS units the orientation is irrelevant.



Suitable for use with ear level receivers.



Shoes not suitable for use with ear level receivers.

Compatibility with direct input leads

By each shoe we also indicate if it is compatible for use with direct input leads.



Suitable for use with direct input leads.




Not suitable for use with direct input leads.

April 2008

OTICON 300P, DigiFocus II Super Power & Spirit 700 SuperPower

Ear level receiver shoe AP201



ID Marks: AP201 Oticon
Contacts: 2 fine contact pins set approx. 1mm apart, beside a black plastic lug
Colour: beige
Manfr. no: 399-50-360-05
Comment: there is a small volume control on one side which should be set to no. 3 when used with CRM220 and no. 4 with fmGenie equipment.

Connevens Pt.No. DOAP201
Price each: £ 17.55

OTICON Personic, Primofocus, MultiFocus Compact, DigiLife and DigiFocus

Audio i/p shoe AP401



ID Marks: AP401 Oticon
Contacts: 2
Colour: translucent
Manfr. no: 399-50-276-02
Comment: there is a small volume control which should be set to no. 3 when used with CRM220 and no. 4 with fmGenie equipment.

Connevens Pt.No. DOAP401
Price each: £ 17.55

Ear level receiver shoe FM4



ID Marks: Compact FM4 Oticon
Contacts: 2 pins, exactly the same as the DAP401, plus a sprung battery contact pin requiring modification of battery door using tool supplied.
Colour: translucent
Manfr. no: 399-51-130-08
Comment: this shoe fits all Oticon compact aids for use with Phonak MLxS units, with the exception of Primofocus Pro2 which is not MLxS compatible.

Connevens Pt.No. DOFM4
Price each: £ 18.20
Supplied with battery door hole tool

OTICON DigiFocus II V2, Adapto, Atlas Plus Direct, GO Direct, Gaia & Spirit II Direct

Audio i/p shoe AP501



ID Marks: AP501 DAI V2 Oticon
Contacts: 1 side contact and 1 sprung pin in base
Colour: translucent
Manfr. no: 399-51-440-03
Comment: there is a small volume control which should be set to no. 3 when used with CRM220 and no. 4 with fmGenie equipment.

Connevens Pt.No. DOAP501
Price each: £ 21.45
Supplied with battery door hole tool

Ear level receiver shoe FM5



ID Marks: Compact FM5 V2 Oticon
Contacts: 2 side contact and 1 sprung pin in base
Colour: translucent
Manfr. no: 399-51-170-07
Comment: Non suffix V2 DigiLife.com & DigiFocus II require factory retro-modification for use with MLx.

MLx: for Oticon WTC compact aids
Connevens Pt.No. DOFM5
Price each: £ 24.69
Supplied with battery door hole tool

OTICON 300P, 380P, 390PL, DigiFocus II SuperPower & Spirit 700 SuperPower

Ear level receiver shoe FM3



ID Marks: WTC Classic Oticon or Oticon FM3
Contacts: 2 fine contact pins set approx. 2mm apart beside a black plastic lug, plus a sprung battery contact pin requiring modification of battery door using tool supplied.
Colour: translucent
Manfr. no: 399-51-180-09
Comment: this shoe fits all current Oticon classic aids for use with Phonak MLxS units.

Connevens Pt.No. DOFM3
Price each: £ 18.20
Supplied with battery door hole tool

OTICON Atlas, Atlas Plus, GO (not Direct), Ergo, Swift, Spirit II, Spirit II VC, Spirit II Power, Spirit VC

Audio i/p shoe AP600



ID Marks: AP-ED01 Oticon
Contacts: a fine contact array
Colour: translucent with grey base
Manfr. no: 399-50-460-06
Comment: there is a small volume control which should be set to no. 3 when used with CRM220 and no. 4 with fmGenie equipment. Also has a HI-LO control which should be set to LO.

Connevens Pt.No. DOAP600
Price each: £ 27.30
Supplied with adaptor tool

Ear level receiver shoe FM6



ID Marks: Oticon FM6
Contacts: a fine contact array
Colour: translucent with grey base
Manfr. no: 399-50-450-04
Comment: Although the DOFM6 shoe version looks identical to the DOFM8, there are internal wiring differences. Does NOT fit Spirit 700 Superpower

Connevens Pt.No. DOFM6
Price each: £ 32.50
Supplied with adaptor tool

OTICON Sumo DM, Sumo XP, Sumo E, Spirit 3 Superpower, Swift 120+

Audio i/p shoe AP700



ID Marks: Oticon DAI AP-700
Contacts: 3 contacts on top
Colour: grey
Manfr. no: 399-50-900-08
Comment: there is a small volume control which should be set to no. 3 when used with CRM220 and no. 4 with fmGenie equipment. Also has a HI-LO control which should be set to LO.

This shoe attaches to the hearing aid by sliding onto the battery door.
Connevens Pt.No. DOAP700
Price each: £ 17.55


Ear level receiver shoe FM7



ID Marks: Oticon FM7
Contacts: 3 contacts on top
Colour: grey
Manfr. no: 399-50-920-02
Comment: This shoe attaches to the hearing aid by sliding onto the battery door.

Connevens Pt.No. DOFM7
Price each: £ 17.55

Digital Hearing Aids for Beginners – a powerpoint presentation on Connevens CD, free if requested with order.



OTICON Spirit 3, Spirit 3 Direct, Spirit 3 VC & Spirit 3 Power, Syncro, Safran, Tego Pro & Tego

Audio i/p shoe AP800



ID Marks: Oticon AP800
Contacts: a fine contact array
Colour: translucent with grey base
Manfr. no: 399-51-480-00
Comment: there is a small volume control which should be set to no. 3 when used with CRM220 and no. 4 with fmGenie equipment. Also has a HI-LO control which should be set to LO.

Connevens Pt.No. DOAP800
Price each: £ 21.97
Does NOT fit Spirit 3 Superpower which requires DOAP700 shoe.
Supplied with adaptor tool

Ear level receiver shoe FM8



ID Marks: Oticon FM8
Contacts: a fine contact array
Colour: translucent with grey base
Manfr. no: 399-51-470-00

Connevens Pt.No. DOFM8
Price each: £ 18.20
Although the DOFM8 shoe version looks identical to the DOFM6, there are internal wiring differences. Does NOT fit Spirit 3 Superpower which requires DOFM7 shoe.
Supplied with adaptor tool

PHONAK Superfront

Audio i/p shoe AS3

BOTH MODELS
Contacts: 3 on one side.
Comment: supplied with a small 3 legged contact plate which may need to be fitted to the hearing aid



Standard Hearing aid microphone live
ID Marks: AS3
Colour: clear with very dark blue base
Manfr. no: 054-0066
Connevens Pt.No. (Standard)
Price each: £ 9.69

Muted Hearing aid microphone muted
ID Marks: MO3
Colour: clear with black base
Manfr. no: 054-0096
Connevens Pt.No. (Muted)
Price each: £ 9.69

Ear level receiver shoe AS3-MLx



ID Marks: AS3-MLx
Colour: clear with a light grey base
Manfr. no: 052-0305
Contacts: 3 contacts on one side, plus an additional battery contact pin for use with MLxS receivers.
Connevens Pt.No. DAS3C
Price each: £ 11.70

PHONAK Audinet, Varionet, Classica & Supero

Audio i/p shoe AS4

BOTH MODELS
Contacts: 3 on one side
Comment: supplied with a small 3 legged contact plate which may need to be fitted to the hearing aid.



Standard Hearing aid microphone live
ID Marks: AS4
Colour: clear with very dark blue base
Manfr. no: 054-0067
Connevens Pt.No. (Standard)
Price each: £ 9.69

Muted Hearing aid microphone muted
ID Marks: MO4
Colour: clear with black base
Manfr. no: 054-0097
Connevens Pt.No. (Muted)
Price each: £ 9.69

Ear level receiver shoe AS4-MLx



ID Marks: AS4-MLxS
Colour: clear with a light grey base
Manfr. no: 052-0306
Contacts: 3 contacts on one side, plus an additional battery contact pin for use with MLxS receivers.
Connevens Pt.No. DAS4M
Price each: £ 11.70

PHONAK Pico, Pico-Forte, Aero, PowerMAXX & Solo T+

Audio i/p shoe AS5A

BOTH MODELS Contacts: 3 on one side



Standard Hearing aid microphone live
ID Marks: AS5A
Colour: clear with very dark blue base
Manfr. no: 054-0092
Connevens Pt.No. DAS5AC
Price each: £ 9.69

Muted Hearing aid microphone muted
ID Marks: MO5A
Colour: clear with black base
Manfr. no: 054-0098
Connevens Pt.No. DAS5AX
Price each: £ 9.69

Ear level receiver shoe AS5-MLx



ID Marks: AS5-MLx
Colour: clear with a light grey base
Manfr. no: 052-0307
Contacts: 3 contacts on one side, plus an additional battery contact pin for use with MLxS receivers.
Connevens Pt.No. DAS5M
Price each: £ 11.70

PHONAK Savia, Eterna, eXtra, Eleva & Una

Audio i/p shoe AS9



General d/i use & MLxS
ID Marks: AS9-MLxS
Contacts: 3 contacts
Colour: beige or palladium (silver)
Manfr. no: 052-3363-0001 (beige)
052-3363-00040 (palladium)
Connevens Pt.No. Beige: DAS9FMB
Palladium: DAS9FMP
Price each: £ 12.80
Supplied with 2 battery drawers and fitting tool

PHONAK Piconet & Sono-Forte P3 (PICS series)

Audio i/p shoe AS6A



ID Marks: AS6
Contacts: 3 on one side
Colour: clear with dark blue base
Manfr. no: 054-0111
Apart from their markings, the difference from the AS5A shoe is in the internal wiring.
Connevens Pt.No. DAS6C
Price each: £ 9.69

Ear level receiver shoe AS6-MLx



ID Marks: AS6-MLx
Colour: clear with a light grey base
Manfr. no: 052-0308
Contacts: 3 contacts on one side, plus an additional battery contact pin for use with MLxS receivers.
Connevens Pt.No. DAS6M
Price each: £ 11.70

PHONAK Nova Forte & P4 Power Zoom

Audio i/p shoe AS7



ID Marks: AS7
Contacts: 4 contacts, 2 either end of the shoe.
Colour: clear with dark blue base
Manfr. no: 054-0113
Connevens Pt.No. DAS7C
Price each: £ 9.69

Ear level receiver shoe AS7-MLx



ID Marks: AS7-MLx
Contacts: 4 contacts, 2 either end of the shoe.
Colour: clear with light grey base
Manfr. no: 052-0309
Although the FM shoe version looks identical there are internal wiring differences.
Connevens Pt.No. DAS7M
Price each: £ 11.70

PHONAK Naída & Exélia

Audio i/p shoe AS10



General d/i use & MLxS
ID Marks: none
Contacts: 3 contacts
Colour: beige, black or taupe
Manfr. no: 052-0022-01 (beige)
052-0022-06 (black)
052-0022-35 (taupe)
Connevens Pt.No. Beige: DAS10BE
Black: DAS10BK
Taupe: DAS10TA
Price each: £ 14.14
Supplied with a replacement battery door unit (shown in taupe) which locks the shoe securely in place.

Connevens Guide to Hearing Aid Direct Input Shoes – 2

Some hearing aids have a direct connection facility, usually referred to as direct input. This facility allows the direct connection of the hearing aid to other audio equipment via a connector usually referred to as a shoe; thus the name 'Hearing aid direct input shoe'. Connevens is able to provide, usually from stock, a wide range of hearing aid direct input shoes; many customers find the number of options overwhelming!



You can put your shoe against the page and identify most shoes without difficulty, making re-ordering easier for you and us!

If you require further help, please do phone our radio aid helpline for assistance when ordering. Please note that although we have usually shown two shoes in the picture to give different views, the shoes are all sold as singles.




Information on compatibility with ear level receivers

By each shoe we indicate whether it is compatible for use with 'ear level receivers' such as MicroMLxS, MLx, MLxS, Lexis, Amigo, Scola and fmGenie WDI units. For use with an MLxS, we also show whether the MLxS needs to be set to a crosswise or lengthwise orientation. For use with fmGenie WDI and MicroMLxS units the orientation is irrelevant.



Compatibility with direct input leads


By each shoe we also indicate if it is compatible for use with direct input leads.



April 2008

SIEMENS Select & Prisma 2M

Audio i/p shoe




ID marks: Labelled A1 inside shoe
Contacts: 3 long fingers
Colour: beige
Manfr. no: 7163178
Comment: distinctive angular shape

Note: Factory retro-modification of the hearing aid may be required for use with MLxS receivers.
Connevens Pt.No. **DAMSELECT**
Price each: **£ 13.00**

SIEMENS Signia & Prisma 2VC+

Audio i/p shoe




ID Marks: none on the outside, but 'for FM' or 'F3' is marked inside in the bottom.
Contacts: four
Colour: translucent
Manfr. no: 5630350

Connevens Pt.No. **DSSIG**
Price each: **£ 16.09**

SIEMENS Swing S1, AMD12+, Selectra, Prisma 2, 2P & Power

Audio i/p shoe



ID marks: C1 inside shoe
Contacts: 3 long fingers
Colour: beige
Manfr. no: 2943970
Note: Factory retro-modification of the hearing aid is required for use with MLxS receivers.
Helpful note: some shoes have a slot in them.

Connevens Pt.No. **DASS1**
Price each: **£ 10.73**

SIEMENS Swing S3 & Prisma 2SP

Audio i/p shoe



ID marks: E1 inside shoe
Contacts: 3 long fingers
Colour: beige
Manfr. no: 2940307

Connevens Pt.No. **DASS3**
Price each: **£ 13.00**

SIEMENS Triano SP, Infiniti Pro SP, Music Pro SP, Prisma 2DSP & Prisma 2DSP+

Audio i/p shoe



ID marks: T-SP or H5 inside
Contacts: 4 contact pins
Colour: translucent
Manfr. no: 7213296

Connevens Pt.No. **DASTRI**
Price each: **£ 13.00**


Hearing aids can be used for pleasure too!

You can use a direct input shoe and a personal stereo direct input lead anywhere you see a headphone socket to link into your hearing aids – see 2008 Connevens Catalogue page 189.



SIEMENS Prisma 2M+

Audio i/p shoe




ID Marks: Prisma 2M+ outside shoe
A5 inside shoe
Contacts: 3 long fingers
Colour: beige
Manfr. no: 8775590
Comment: distinctive angular shape, physically the same as the DAMSELECT shoe but not compatible.
Note: When fitting this shoe a firm push to click fully home is required.

Connevens Pt.No. **DSIEA**
Price each: **£ 13.00**

SIEMENS Infiniti Pro

Audio i/p shoe



ID Marks: G5 inside shoe
Contacts: a fine array of 4 contacts
Colour: translucent
Manfr. no: 7770584
Comment: this shoe is identical to DSIEC (you will notice they are both ID'd G5), except for the colour coded contact covers for the hearing aid which are model specific.

Connevens Pt.No. **DSIEB**
Price each: **£ 13.00**

SIEMENS Prisma Pro & 2 Pro

Audio i/p shoe




ID marks: G5 inside shoe
Contacts: a fine array of 4 contacts
Colour: translucent
Manfr. no: 10056036
Comment: this shoe is identical to DSIEB (you will notice they are both ID'd G5), except for the colour coded contact covers for the hearing aid which are model specific.

Connevens Pt.No. **DSIEC**
Price each: **£ 13.00**

SIEMENS Acuris S & P, Cielo 2, Artis 2S, Artis 2SVC, Centra P, S & SVC, Triano S & SL

Audio i/p shoe




ID marks: None
Contacts: four
Colour: translucent
Manfr. no: 7220424

Connevens Pt.No. **DSIED**
Price each: **£ 20.80**

SIEMENS Chroma VC & P, Centra HP, Artis S, P & 2P, Cielo P & 2P, Prisma 2 Pro P & DP & Reflex DP

Audio i/p shoe




ID Marks: L5 inside shoe
Contacts: 4 on one side
Colour: translucent
Manfr. no: 10054830

Connevens Pt.No. **DSIEG**
Price each: **£ 13.00**

SIEMENS Prisma 2 Pro DSP & SP, Artis SP & 2SP, Centra SP, Cielo 2SP, Intuis SP DIR & Chroma SP

Audio i/p shoe




ID marks: M5 inside shoe
Contacts: 4 on one side
Colour: translucent
Manfr. no: 10054617

Connevens Pt.No. **DSIEI**
Price each: **£ 13.00**

GN ReSound Canta 780D, 280D, Danalogic 283D

Audio i/p shoe DA11




ID Marks: GN ReSound DAI 11
Contacts: 4 in two pairs
Colour: translucent
Manfr. no: 15072200

Connevens Pt.No. **DGNR780D**
Price each: **£ 10.14**

GN ReSound Danalogic6 (6060 & 4060), Azure (AZ60), Metrix (MX60), Pixel (PL60) & Plus5 (RP60)

Audio i/p shoe



ID Marks: none
Contacts: 4 finger contacts on one side
Colour: translucent
Manfr. no: 15832600

Connevens Pt.No. **DGNRA**
Price each: **£ 10.14**

Even more shoes available online...

In addition to the shoes shown here, Connevens also has stocks of many shoes for non current hearing aids and the website is regularly updated with details of newly available shoes.


www.DeafEquipment.co.uk

Use the search engine to find the correct shoe for your aid – fill in the make or model of your hearing aid.



UNITRON ICON, Sound FX & BE39FX series

Audio i/p shoe




ID Marks: none, but there is a Unित्रon trademark.
Contacts: 4
Colour: translucent
Manfr. no: 054-5031-13
Comment: there is a microphone mute switch on the edge of the shoe.

Connevens Pt.No. **DUI**
Price each: **£ 13.81**

UNITRON US series

Audio i/p shoe




ID Marks: Unित्रon trademark.
Contacts: 4
Colour: translucent
Manfr. no: 054-5026-13
Comment: there is a microphone mute switch on the edge of the shoe.

Connevens Pt.No. **DUUS**
Price each: **£ 16.25**

UNITRON Unison 3HP, Unison 6HP & Essential HP

Audio i/p shoe



ID Marks: Unित्रon hearing 384-0A
Contacts: 2 contact, 1 pin
Colour: translucent
Manfr. no: 054-5119-13
Comment: This is a combined audio i/p shoe and battery tray. The user requires average DIY as this 'shoe' replaces the existing battery drawer with this battery tray which incorporates the sockets for the d/i or ear level receiver.

Connevens Pt.No. **DUAIBT**
Price each: **£ 13.81**

UNITRON Unison 3 & 3P, 6 & 6P, Essential Standard & Essential P

Audio i/p shoe



ID Marks: Unित्रon hearing 384-0A
Contacts: 2 contact, 1 pin
Colour: translucent
Manfr. no: 054-5123-13
Comment: the user requires average DIY as this 'shoe' replaces the existing battery drawer with this battery tray which incorporates the sockets for the d/i or ear level receiver.

Connevens Pt.No. **DUAIBT13**
Price each: **£ 13.81**

UNITRON 2, 4, Unison 2P & Conversa

Audio i/p shoe ASSA



These Unित्रon Hearing aids use the three variations (same MLxS compatibility etc.) of Phonak ASSA shoe – please see Connevens Catalogue 2008 page 69 for pricing and part numbers.

Visit the website!

New direct input shoes updated regularly plus some older shoes are still available for online ordering.


We even have video clips showing what to do with some of the extra bits and bobs and tools!



Go to www.DeafEquipment.co.uk select Hearing Aid Direct Input Shoes and follow the link to watch video clips.

WIDEX Senso Diva SD-19, Inteo IN-19 and Vita SV-19

Audio i/p shoe Widex 5




ID Marks: Widex 5
Contacts: 5
Colour: translucent
Manfr. no: 295 W5

Connevens Pt.No. **DW5**
Price each: **£ 14.63**

WIDEX for Senso Vita SV-38, P37/P38, Bravissimo BV-38, Bravo B32 and Logo L32

Audio i/p shoe Widex 6



ID Marks: Widex 6 (previously Widex 4, 3i, 3 and before that Logo Senso FM)
Contacts: 3 on one side, 1 on the other
Colour: translucent
Manfr. no: 295 DW6
Comment: Has small catch on one side to lock shoe in place. (size 675 battery hearing aid) With a slide on/off connection for the Widex SCOLA fm receiver

Connevens Pt.No. **DW6**
Price each: **£ 14.63**

WIDEX Senso Diva SD9, Senso P7/P8, Inteo IN-9E Bravo B1/B2 & Vita SV-9

Audio i/p shoe Widex 7



ID Marks: Widex 7 (previously Widex 1)
Contacts: 4, the centre left pin is offset and angled
Colour: translucent
Manfr. no: IN9
Comment: Size 13 battery hearing aid. With a slide on/off connection for the Widex SCOLA fm receiver

Connevens Pt.No. **DW7**
Price each: **£ 14.63**

WIDEX Logo Senso C series and Bravo B11/B12 and Logo L6-L12

Audio i/p shoe Widex 8



ID Marks: Widex 8
Contacts: 3 on one side, 1 on the other
Colour: translucent
Manfr. no: 295 W8
Comment: Fits Senso C8+ - C19+. Has small catch on one side to lock shoe in place. (size 13 battery hearing aids). With a slide on/off connection for the Widex SCOLA fm receiver

Connevens Pt. No. **DW8**
Price each: **£ 14.63**



Connevens Limited
EQUIPMENT FOR DEAF AND HARD OF HEARING PEOPLE
IN EDUCATION, EMPLOYMENT, THE HOME & LEISURE

54 Albert Road North, Reigate, Surrey RH2 9YR, United Kingdom
Switchboard: 01737 247571 Minicom: 01737 243134 Fax: 01737 223475

Information website www.connevens.com
Online shopping website www.DeafEquipment.co.uk