



1a

fmGenie



# WDI – Wireless Direct Input fmGenie system – QFMG8

User's own hearing aids and direct input shoes

## For individuals using direct input hearing aids.

The QFMG8 package includes the fmGenie equipment, accessories and recommended spares – all you need to add are the appropriate connecting shoes – for a working direct input system.



### QFMG8 system package includes:



fmGenie transmitter	Part no. FMG110	£ 459.35
Microphone pack (600mm lead)	FMG213M6	46.44
Spare 600mm microphone lead	FMG21M600	10.00
fmGenie receiver	FMG121	525.12
WDI unit (qty two)	FMGWDI	106.04
Screwdriver set	*Y015*	2.00
500mm WDI neck loop (qty two – inc. one as a spare)	FMG474M	53.04
Radio aid carry pouch	FMG615	7.95
Pouch with belt clip (qty two)	FMG630	21.72
Hip pouch & shoulder strap	FMG640	12.48
Transmitter neck harness & pouch	FMG621	12.48
AA Hybrio ready use rechargeable batteries (pk of 4)	MBRHAAX4-	4.61
Twin tail charger	FMG520	19.62
Listening test set	FMG711	8.34
fmGenie operating instructions – includes fast track guide and laminated check sheets	–	–
Shipping and order processing	PSH-PARCEL	9.95

Please note, this package does not include direct input shoes – these must be purchased separately, see Section 2.

<b>QFMG8 Direct input fmGenie system</b>	<b>£ 1299.14</b>
<b>QFMG8 with 2 year warranty</b> – QFMG8 with additional year 2 warranty	<b>£ 1398.14</b>
<b>QFMG8 with 3 year warranty</b> – QFMG8 with additional years 2 & 3 warranty	<b>£ 1497.14</b>

## A few fmGenie Questions & Answers

**I have an fmGenie direct input package (QFMG3), I find it very good in one-to-one situations but not so good in group discussions, any suggestions?**

Yes, add an fmGenie table top conference microphone (FMG220).

**Why do you not suggest using direct input for the QFMG5 fmGenie visitor's system?**

As different hearing aids require different shoes you would be unlikely to have the correct ones!

**Will a hearing aid user hear any sound quality difference between conventional direct input and WDI?**

No.

**Do I use my hearing aid on 'T' when using the WDI neck loop?**

No, the fmGenie WDI facility has nothing to do with hearing aid 'T' pickup.

**I know that my hearing aid can be used with MLxS units is it also suitable for use with an fmGenie WDI unit?**

Yes.

**I have some spare fm (MLx) compatible shoes – are they suitable for use with the WDI?**

Yes.

**My fmGenie system is quite old and I want to use it with**

**the WDI, do I need to send it back for modification?**

No, using a WDI unit is easy. All fmGenie receivers, past or present, can be used with WDI – all you need to do is to change the settings in the receiver as explained in the instructions.

**Can I use my present 'T' neck loop with a WDI?**

No, you must use the special WDI neck loop which is different.

**The WDI neck loop has the same plug that is used on my Walkman, will it work with it?**

No. The fmGenie receiver has special circuitry in it to drive the WDI neck loop. Although the WDI is a recent product addition the WDI drive circuitry has always been in the fmGenie receiver. When you change the fmGenie receiver settings and add a link plug you activate this special circuitry.

**Can I use an fmGenie tx with ear level fm receivers?**

Yes, you need to turn the transmitter range to LR (Long Range).

*Many more FAQ's can be found in the fmGenie User Guide & Reference Manual*

