

## Baha spares and accessories

In conjunction with Cochlear Europe, Connevens are delighted to offer a Baha spares and accessories supply source.

We have a direct link web address where you can find more information and order online – [www.bit4bahas.co.uk](http://www.bit4bahas.co.uk).



Also please visit here to view our full range of products – if the item you require is not there please email [sales@connevens.com](mailto:sales@connevens.com) and we will try to source it for you.

### Telecoil unit for Baha Divino and Compact

Allows Baha users to enjoy the benefits of inductive listening. The telecoil unit has a 3 setting switch: 'M' microphone, 'M+T' microphone & telecoil mixed and 'T' telecoil.

Size: 56 x 8mm Weight: 2g.

Manfr. no: 90185

DCBA185


£ 180.60 each



### Snap coupling abutment covers



The abutment cover can be attached to the abutment when the Baha sound processor is not in place. It gives an aesthetically pleasing appearance to the abutment and protects the inside of the abutment from dirt. The cover is supplied with the sound processor kits at the time of implant, but we offer it for purchase separately.



### For Compact/Classic 300/Cordell II £ 9.00 each

Manfr. no: 90751 (Beige), 90752 (Black), 90753 (Grey)

DCBA751 Beige

DCBA752 Black

DCBA753 Grey

### For Divino £ 9.00 each

Manfr. no: 90744 (Blonde), 90745 (Black), 90746 (Silver), 90747 (Brown)

DCBA744 Blonde

DCBA745 Black

DCBA746 Silver grey

DCBA747 Brown

### Phonak MLxS Baha receiver module

MLxS units for use with a Baha have a different connection pin configuration. The MLxS receiver connects directly into the Baha unit. A MicroMLxS version is not available for Bahas.

3PMLXSBA

£ 340.00 each



### Baha MLxS tester

Adaptor to allow testing of a Baha MLxS with an FMG750 Accessories Tester or a direct i/p hearing aid.

3PMLXSBTES

£ 75.00 each




### Baha soft cleaning brush

A replacement soft cleaning brush for use with Bahas – it is important that you keep your Baha and the skin around the abutment clean to prevent soreness and irritation.

Manfr. no: HCB738

DCBA738

£ 2.15 each



### Caring for the abutment area

**Daily care:** Clean the skin around the abutment daily to remove any skin debris. A soft cleaning brush is given when the sound processor is fitted. Use it with soap and water to keep the area clean.

**Weekly care:** Skin debris may also build up inside the abutment or on the abutment snap itself. This should also be removed by using the soft cleaning brush.

**Regular inspection of the abutment:** always check that the abutment is firmly in place. This can be done yourself or ask your clinic to help you during your routine appointments, about every six months.

*Contact your doctor if you have persistent problems such as soreness, infection, or inflammation around the abutment.*

### Baha safety line

Discreet clothes safety line to help prevent loss or damage of your Baha.

Manfr. no: 90711

DCBA711

£4.00 each



### Baha softbands

Intended for babies and small children who can benefit from a Baha sound processor but who are, as yet, too young for surgery. We supply these Softbands as spares or replacements for those issued by your hearing professional.



The softband combines an adjustable elastic band with a snap connector disk which is pressed against the skin behind the ear, or at another bony location of the skull, through the pressure from the band. The sound processor is attached to the disk and the sound is then transmitted through the bone of the skull to the functioning inner ears.

*The softband comes in a kit containing two softbands.*

Manfr. no: 90544 / 90545 / 90851



DCBA544 1 x space design & 1 x plain blue £ 70.50 each

DCBA545 1 x teddy design & 1 x plain red £ 70.50 each

DCBA851 1 x pale pink & 1 x pale blue £ 70.50 each