

# Hearing Aid Direct Input Shoes



2

Hearing Aid Direct Input Shoes

Some hearing aids have a direct connection facility, usually referred to as direct input. This facility allows the direct connection of the hearing aid to other audio equipment via a connector usually referred to as a shoe; thus the name 'Hearing aid direct input shoe'.

Connevens is able to provide, usually from stock, a wide range of hearing aid direct input shoes; many customers find the number of options overwhelming! You can put your shoe against the page and identify most shoes without difficulty, making re-ordering easier for you and us!

If you require further help, please do phone our radio aid helpline for assistance when ordering. Please note that although we have usually shown two shoes in the picture to give different views, the shoes are all sold as singles.



**VAT** All direct input shoes: VAT relief is permitted subject to the receipt of a valid VAT declaration form.

## Information on compatibility with ear level receivers



By each shoe we indicate whether it is compatible for use with 'ear level receivers' such as MicroMLxS, MLx, MLxS, Lexis, Amigo, Scola and fmGenie WDI units.

For use with an MLxS, we also show whether the MLxS needs to be set to a crosswise or lengthwise orientation. For use with fmGenie WDI and MicroMLxS units the orientation is irrelevant.



Suitable for use with ear level receivers.



Shoes not suitable for use with ear level receivers.

## Compatibility with direct input leads

By each shoe we also indicate if it is compatible for use with direct input leads.




Suitable for use with direct input leads.



Not suitable for use with direct input leads.

## GN ReSound Canta 780D, 280D, Danalogic 283D

Audio i/p shoe DA11



ID Marks: GN ReSound DAI 11  
Contacts: 4 in two pairs  
Colour: translucent  
Manfr. no: 15072200  
Connevens Pt.No.  
Price each:

DGNR780D  
£ 10.14

## GN ReSound Danalogic6 (6060 & 4060), Azure (AZ60), Metrix (MX60), Pixel (PL60) & Plus5 (RP60)

Audio i/p shoe



ID Marks: none  
Contacts: 4 finger contacts on one side  
Colour: translucent  
Manfr. no: 15832600  
Connevens Pt.No.  
Price each:

DGNRA  
£ 10.14

## Visit the website!

New direct input shoes updated regularly plus some older shoes are still available for online ordering. We even have video clips showing what to do with some of the extra bits and bobs and tools!



Go to [www.DeafEquipment.co.uk](http://www.DeafEquipment.co.uk) select **Hearing Aid Direct Input Shoes** and follow the link to watch video clips.

## Hearing aids can be used for pleasure too!

You can use a direct input shoe and a personal stereo direct input lead anywhere you see a headphone socket to link into your hearing aids – see page 189.

